

National Freezer Challenge

1 September - 31 December 2025

We encourage laboratories to clean out their freezers (-20°C, -70/-80°C and cryogenic freezers) and use cold storage space more efficiently.

One single ULT freezer uses as much electricity as 2-3 average Dutch households!


Together, we can create a more sustainable lab environment and save energy.

Participate? Mail to
[mailadres]
or scan the QR code below!


Freezer
-20°C




Count the amount of cleared boxes




ULT freezer
-70°C or -80°C



Count the amount of cleared boxes



Cryogenic Freezers
-150°C and Liquid Nitrogen



Count the amount of cleared samples



More info

